

**Amendments to the Specification:**

Please replace the paragraph beginning on Page 6, Line 40 and ending on line 44 with the following amended version of this paragraph:

Figs. 20 and 21 are pictorial and end views, respectively, of a forearm supination and pronation device 1300 which includes a base 1301, right and left adjustable base flanges 1302L, 1302R, respectively, an inflatable member 1320, straps 1325, a retaining cuff 1327, and a glove assembly [[1350]] 1340 including a thumb retainer 1351.

Please replace the paragraph beginning on Page 14, Line 31 and ending on line 36 with the following amended version of this paragraph:

This apparatus 900 includes a substantially flat base 901, an inflatable member 920, and a floating yoke member 930. The substantially flat base 901 includes a pair of upwardly directed holes 902 (as Fig. 14 is viewed) which slidably accept a corresponding pair of downwardly-directed feet [[911]] 931 which extend downwardly and leftwardly at an incline from the inclined yoke plate 910.

Please replace the four paragraphs beginning on Page 17, Line 29 and ending on Page 18 Line 6 with the following amended version of these paragraphs:

Figs. 20 and 21 are pictorial and end views, respectively, of a forearm supination and pronation device 1300. This device 1300 includes a base 1301, right and left adjustable base flanges 1302L, 1302R, respectively, an inflatable member 1320, straps 1325, a retaining cuff 1327, and a glove assembly [[1350]] 1340 including a thumb retainer 1351. The base 1302 is configured to be stationary and to provide support for the other members.

The right and left adjustable base flanges 1302L, 1302R, respectively, are configured to be attached to the base, and to be adjustable. During use they are rigidly affixed relative to the base [[1302]] 1301.

The straps 1325 are positioned intermediate one of the base flanges and the rigid post member 1350.

Attorney Dkt: 45011-271086

The rigid post member 1350 is attached to the glove assembly 1340. This attachment can be temporary through VELCRO, snaps, or even a clothespin configuration. Alternately the post 1350 could be pivotally attached to the **[[plate]]** base 1301.

Attorney Dkt: 45011-271086